# User manual appendix B

**myopro**<sup>®</sup>  
motion w

*Specific instructions for your unique MyoPro 2*

**myomo**<sup>®</sup>  
my own motion


## CONTENTS

› MyoPro 2 Motion W Diagram	5
› Putting Your MyoPro 2 On (“Donning”)	7
› Steps for Putting Your MyoPro 2 On	10
› Taking Your MyoPro 2 Off (“Doffing”)	11
› Steps for Taking Your MyoPro 2 Off	13
› Using Your MyoPro 2	15
› Control panel	17
› Modes of Use	19
› Start Mode	19
› Elbow Modes	19
› Changing elbow Modes with the Control Panel	22
Dansk	23
Deutsch	46

## MYOPRO 2 MOTION W DIAGRAM

Component names and references (not every device includes every component).

# PUTTING YOUR MYOPRO 2 ON ("DONNING")

## PUTTING YOUR MYOPRO 2 ON (“DONNING”)

The following pages include step-by-step instructions for putting the MyoPro 2 on (sometimes referred to as “donning the MyoPro 2”).

### **BEFORE donning the MyoPro 2, review the warnings below:**

-  **WARNING:** DO NOT expose the MyoPro 2 to flame or excessive heat; personal injury may occur.
-  **WARNING:** The MyoPro 2 is a prescription device and should only be used as indicated under clinical supervision of a Provider, or after receiving instruction from a clinical Provider for at-home use.
-  **WARNING:** The MyoPro 2 is only to be used by the person for whom it is prescribed.
-  **WARNING:** Each component of the device is supplied for safety and best performance. Use all provided components for optimal operation.
-  **WARNING:** ONLY use the MyoPro 2 upper limb orthosis on the prescribed elbow, wrist and hand joints. Each MyoPro 2 is built specifically for a left arm or a right arm; never wear the MyoPro 2 on the other arm.
-  **WARNING:** Tight straps may restrict the User’s circulation. Therefore, always check that straps are not too tight throughout the User’s range of motion.
-  **WARNING:** Do not sleep while wearing the device as straps may alter the User’s circulation.
-  **CAUTION:** If you detect fumes, flames, melting of components around the battery or battery charger, or the device is hot to touch, TURN THE DEVICE OFF IMMEDIATELY and contact your Provider.

## STEPS FOR PUTTING YOUR MYOPRO 2 ON

1. **POWER OFF:** While donning, the MyoPro 2 should be turned OFF; check that it is OFF before starting to put the MyoPro 2 on your arm.
  - ⌚ **NOTE:** Pads and straps should be fully air-dried from the previous use.
  - ▲ **CAUTION:** Do not use a hairdryer to dry components of the MyoPro 2.
2. **PREPARE:** Place the MyoPro 2 in your lap, or on a flat surface that is not slippery.
3. **HAND:** Start by feeding the fingers and hand through the Hand Support Shell. Be sure the strap around the support shell sits at the base of your thumb. Tighten the Velcro strap around the hand.
  - ⌚ **NOTE:** The motor and battery case should rest on the outside of the arm. If the motor and battery case is aligned along the interior of the arm, you may be donning upside down or on the wrong arm.
4. **WRIST:** Lower your wrist into Multi-Articulating Wrist (MAW), Rotation Portion
5. **FOREARM:** Lower your forearm into the Forearm Shell. Secure the Velcro of the Forearm Closure.
6. **UPPER ARM:**
  - Wrap the Upper Sensor Cuff around your upper arm, keeping attention on the Biceps and Triceps sensors; when the Upper Sensor Cuff is properly fit, the sensors should be centered over the biceps and triceps muscles.
  - ⌚ **NOTE:** the sensor with the **BLUE** Velcro goes on biceps, the sensor with the **RED** Velcro goes on triceps.

- Secure the Upper Cuff Closure.
7. **HARNESS:** Secure the harness:
    - Rest the shoulder pad on the shoulder of your MyoPro 2 arm.
    - Wrap the strap under your non-MyoPro 2 arm.
    - Secure quick release tab.

# TAKING THE MyoPro 2 OFF ("doffing")

## TAKING THE MYOPRO 2 OFF (“DOFFING”)

 **WARNING:** If at any time during the use of this device, you notices any of the following, discontinue use and seek guidance from your Provider.

- Movement does not match the User’s desired motion.
- Persistent redness, swelling, or skin breakdown (bleeding, chafing, etc.)
- Rash on the arm, hand, or fingers.
- Pain associated with wearing the MyoPro 2 orthosis.
- Unusual noises from the orthosis (popping, clicking, etc.)
- Smells from the orthosis (smoking, burning plastic, etc.)
- Odor from the orthosis (sour smells or other indications of bio-contamination.)

## STEPS FOR TAKING YOUR MYOPRO 2 OFF

1. Turn the MyoPro 2 OFF.
2. Place the MyoPro 2 in your lap, or on a flat surface that is not slippery.
3. Undo the harness Quick Release tab.
4. Undo the Upper Cuff Closure.
5. Undo the Forearm Closure.
6. Gently slide your hand back through the loosened hand strap and Hand Support Shell
7. Rest the MyoPro 2 in a clean, dry location to ensure it will fully air-dry before next use.
8. Place the MyoPro 2 into the carrying bag for travel or storage.

 **NOTE:** The MyoPro 2 should remain OFF any time the MyoPro 2 is not being worn.

# USING YOUR MYOPRO 2


## CONTROL PANEL

See the *User Manual – Sections for All MyoPro 2 Devices* document for a comprehensive Control Panel Buttons & Lights Diagram.

⌚ **NOTE: For your MyoPro 2 Motion W**

- The Grasp mode button is deactivated
- The Close Light and Open Light may turn on if a certain error occurs. See MyoPro 2 Error section of the User Manual

## MODES OF USE

 **WARNING:** Use caution when using the device in shoulder positions where it is possible for the User to hit him or herself.

The MyoPro 2 may be used in a variety of **modes**. The different modes allow you to operate the MyoPro 2 by tensing or relaxing different muscles. You may find that you want to use different muscles – and therefore different modes – when performing specific tasks with your MyoPro 2. Modes may be changed using the buttons on the control panel.

When the device is fit to you, your Provider will optimize the sensitivity settings within each mode depending on the strength of your EMG signal at that time. If over time you feel that the MyoPro 2 is not assisting your arm or hand adequately, please contact your Provider; the sensitivity settings may need to be readjusted.

### START MODE

When you first power on the MyoPro 2, the elbow will be in Standby Mode. This setting may be changed by your Provider.

### ELBOW MODES

The elbow has four distinct modes of use when the MyoPro 2 is powered on.

1. Standby mode
2. Biceps mode
3. Triceps mode
4. Dual mode

These are explained in the following pages.

## ELBOW MODES, CONTINUED

<p><b>Standby mode</b></p>	<p>In this mode, neither the Biceps or Triceps Light will be illuminated. The elbow motor will not respond to your EMG signal from either muscle group.</p> <p>Though the MyoPro 2 is powered on, no assistance is being given to your elbow, and the arm will neither flex nor extend. The sensors will be reading your EMG signal from both biceps and triceps, but the motor response is paused.</p> <p>Actions:</p> <ul style="list-style-type: none"> <li>• This can be useful as a resting mode if straps or sensors need to be adjusted, or if you want to take a short break without turning the device off.</li> </ul>	
<p><b>Biceps mode</b></p>	<p>In this mode, the Biceps Light will be illuminated. The elbow motor will respond to your biceps EMG signal.</p> <p>Actions:</p> <ul style="list-style-type: none"> <li>• When you relaxes your biceps, the device will extend.</li> <li>• When you contracts your biceps, the device will flex.</li> </ul>	

## ELBOW MODES, CONTINUED

<p><b>Triceps mode</b></p>	<p>In this mode, the Triceps Light will be illuminated. The elbow motor will respond to your triceps EMG signal.</p> <p>Actions:</p> <ul style="list-style-type: none"> <li>• When you relax your triceps, the elbow will flex.</li> <li>• When you contracts your triceps, the elbow will extend.</li> </ul>	
<p><b>Dual mode</b></p>	<p>In this mode, both the Biceps and Triceps Light will illuminated. The elbow motor will respond to your biceps <i>and</i> triceps EMG signal.</p> <p>The elbow will only respond to muscle contraction (not relaxation, as in other modes) to assist you with active flexion and active extension of your elbow.</p> <p>You must engage both muscle groups in order for your arm to move in both directions. The MyoPro 2 will respond to whichever muscle signal is stronger at a given moment in time. In order for Dual mode to be most effective, you should be able to isolate the desired muscle (biceps or triceps) and have a limited number of co-contractions prior to using Dual mode. When used well, this mode should allow you greater ability to hold an arm position at a mid-point between fully flexed and fully extended.</p> <p>Actions:</p> <ul style="list-style-type: none"> <li>• When you contract your biceps, the elbow will flex, as long as your biceps signal is greater than your triceps signal.</li> <li>• When you contract your triceps, the elbow will extend, as long as your triceps signal is greater than your biceps signal.</li> </ul>	

## CHANGE ELBOW MODE WITH THE CONTROL PANEL

The elbow mode button cycles through the 4 elbow modes, moving to the next mode each time you press the mode button.

**When the MyoPro 2 is turned ON, the elbow will be in Standby mode unless it has been changed by your Provider.**

●	Press the elbow mode button one time to put the elbow into <b>Biceps mode</b>
● ●	Press the elbow mode button again (a second time) to put the elbow into <b>Triceps mode</b>
● ● ●	Press the elbow mode button again (a third time) to put the elbow in <b>Dual mode</b>
● ● ● ●	Press the elbow mode button again (a fourth time) to return the elbow <b>Standby mode</b>

Repeat the above cycle to scroll through and change to another mode at any point during your MyoPro 2 use.

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